

Additional Sizing Information

Front Leg Splints

	X-Large	Large	Medium	Small	X-Small	Feline
	J-634md	J-634mc	J-634mb	J-634ma	J-634m	J-633m
A	12"	11"	10½"	9"	8"	5½"
B	7½"	7"	6"	5¼"	4½"	3"


Rear Leg Splints

	X-Large	Large	Medium	Small	X-Small	Feline
	J-632md	J-632mc	J-632mb	J-632ma	J-632m	J-633ma
A	16"	14"	11"	10"	9½"	6¼"
B	11¼"	9½"	8"	7"	6½"	4¼"

Front Leg Splint

"A" Measurement: Total Length


"B" Measurement: Carpel Pad to the end of Toes



Rear Leg Splint

"A" Measurement: Total Length

"B" Measurement: Hock to Toes



Sizing and Fitting

Feline Size	For cats and small dogs less than 10 lbs. (Dachshund)
X-Small Dog	For dogs 10-25 lbs. (Terriers)
Small Dog	For dogs 25-40 lbs. (Shelties)
Medium Dog	For dogs 40-65 lbs. (Border Collies)
Large Dog	For dogs 65-85 lbs. (Labrador & Golden Retriever)
X-Large Dog	For dogs 85-110 lbs. (Rottweiler)
XX-Large Dog	For dogs 110 lbs. and up (Great Dane) <i>will be available soon</i>