


DIAGNOSTIC ECG and MONITOR

Built in printer



Large easy to read screen


Touch keypads

Easy keypad set up


The ECG is one of the most under utilized diagnostic parameters in veterinary medicine. An ECG is an absolute necessity for proper cardiac diagnostics. The Jorvet Single Channel ECG and monitor offers what veterinarians need at a **very reasonable price**. The all-important lead II rhythm strip can be run for screening purposes; and if further evaluation is deemed necessary, all 6 common leads can be run for more thorough diagnostics.

FEATURES AND BENEFITS

- ♥ One touch keypad buttons allow for easy set up and operation
- ♥ Large screen is well lit and easy to read
- ♥ Single channel displays and prints one lead at a time
- ♥ Features both **MANUAL** and **AUTO** modes
- ♥ Dimensions are 12”L x 8 7/8”W x 2 5/8”H
- ♥ Weighs just 6.25 lbs
- ♥ Operates on either AC or 12 volt NiMH rechargeable battery
- ♥ Patient comfortable electrode clips
- ♥ Built in paper printer
- ♥ Supports *patient monitoring* via onscreen tracing, i.e., lead II
- ♥ Supports *diagnostic evaluation* via the printed tracing
- ♥ **3 year warranty on the ECG unit, 1 year on battery and leads**



Rounded clips are very pet friendly!



includes
BONUS GIFT! Larry Tilley's *Electrocardiography for the Small Animal Practitioner*
 Takes the mystery out of properly using an ECG machine and evaluating the results

J1010—ECG machine

J1010d1—Complete Lead Set

J1010d2—4-Wire Limb Lead Only

J1010d3—ECG paper rolls (each)