



## JorVet Anesthesia Plus

Deluxe version of the standard JorVet Anesthesia machine.

- Plenty of shelves, powdered coated x-ray steel 9" width x 24" long
- Lower shelf 12" x 24" long for storing patient monitors, electrosurgery unit or any other equipment
- Plenty of room for extra vaporizer or additional flow meter
- Handle for stability when moving
- Room for e tank yoke
- Large locking casters for east and smooth motion

Include universal F anesthesia

- 3 liter breathing back
- Non rebreathing anesthesia

**J558N** w/o vaporizer

**J558NA** w/ Tech III refurbished isoflurane or choose from new and refurbished sevoflurane vaporizer



**Jorgensen Laboratories, Inc.**  
1450 Van Buren Ave  
Loveland, CO 80538  
(800) 525-5614  
Info@JorVet.com  
www.JorVet.com