





Clear Adhesive Surgical Drape with Povidine



The majority of human surgeries use clear adhesive drape over the incision site. The rest of the patient is draped with standard non-woven disposable drape material. The clear adhesive drape offers several advantages:

1. **Transparency** - the film is ultra thin and highly transparent, which gives clear observation to the surgical field for locating anatomical landmarks.
2. **Antibacterial** - impregnated with povidine iodine giving the drape a slight **orange** tint. Very effective against most skin pathogens. Emergence of MRSA makes this more imperative.
3. **Fluid and Pathogen resistance** - sterile barrier prevents fluids or pathogens right up to the edge of the incision.
4. **Permeable** - allows the skin to breathe and vapor to pass.
5. **High Elasticity** - follows natural contours for reliable gap free adherence.
6. **Low Sensitivity** - no irritation to skin or wound
7. **Moderate Viscosity** - no injury to skin when removing the dressing.

Adhesive drape at incision's edge.



Slight orange tint.

Range of sizes, packaged sterile in sealed pouch

J1086a	14cm x 20cm (6" x 8")	20 pieces/box
J1086b	20cm x 30cm (8" x 12")	30 pieces/box
J1086c	30cm x 45cm (12" x 18")	20 pieces/box
J1086d	45cm x 60cm (18" x 24")	10 pieces/box

Jorgensen Laboratories, Inc.

1450 Van Buren Avenue Loveland, CO 80538 (800) 525-5614 Info@JorVet.com www.JorVet.com