




Fecal Assay Kits



JorVet in conjunction with Dr. Don Bliss, has developed several species specific fecal assay kits. Dr. Bliss is a noted parasitologist from Wisconsin who has lectured extensively throughout the world on parasite control. These assay kits contain everything a veterinarian or a trained lay person would need to do a large number of fecal exams at a time. It does not include a microscope or centrifuge.

The emerging resistance to some anthelmintic agents requires fecal assays be done. This is an opportunity for veterinarians to provide their professional expertise regarding parasite control.

J1101ED Equine Fecal Assay Kit

- Sheather's sugar flotation solution, 1 gallon
- Dosing syringe for fecal flotation solution
- Microscope slides and coverslips
- Wash bottle
- Test tube rack
- Centrifuge tubes
- Tongue depressors for stirring
- Fecal assay bucket
- Pet mitts for stool collection
- Fecal assay strainer
- Hand held tally counter
- Surgical scissors
- Colored equine parasite poster
- Dixie cups



J1101RD Ruminants Fecal Assay Kit

Same as above only with colored ruminant parasite poster

J1101SD Small Animal Fecal Assay Kit

Same as above only with colored small animal parasite poster
Uses contrate filters instead of a strainer



J-1101ED

Jorgensen Laboratories, Inc.

1450 Van Buren Avenue Loveland, CO 80538 (800) 525-5614 Info@JorVet.com www.JorVet.com