



FitPAWS® Giant Rocker Board

Finally, a Giant Rocker Board designed specifically for dogs! The FitPAWS Giant Rocker Board has an extra-wide training surface 30" x 30" (75cm x 75cm) and a 1.5" (4cm) center fulcrum bar. Depending on how your dog is facing, the board will rock from side-to-side or front-to-front. Use this canine conditioning tool for rehab, balance training, and beginner seesaw and contact training for agility.



- Padded, non-slip textured training surface
- Padding on edges of board to muffle sound when in use on indoor flooring
- Easy-to-clean. Surface mat is attached w/Velcro and is easily removed for cleaning. Hang up mat and spray with water. Allow to dry completely before replacing.
- Easy to maintain. Replacement FitPAWS Mats are available for purchase.
- Includes insert with exercise examples.

Product Code

J1127

J1127m

Description/Size

30" x 30" (75cm x 75cm) Giant Rocker Board

24" x 60" (61cm x 152cm) Replacement Mat

FitPAWS® Hurdle Set



An incredible value for the backyard agility enthusiast, traveling animal rehab practitioner or professional dog trainer. Completely adjustable and portable dog agility cone and hurdle set includes four 40" (102cm) jumping bars, eight highly-visible orange hurdle cones, eight height-adjustment attachments and a black nylon carrying bag.

- Use indoors or outdoors
- Height of jumping bar adjusts in 2" (5cm) increments from 2"-30" (5cm-76cm)
- Perfect for easy rehabilitation or more challenging agility training
- Lightweight for on-the-go portability
- Some assembly required

Product Code

J1126

Description

Hurdle & Cone Set

Lowest setting 2" (5cm)

Highest setting 30" (76cm)





FitPAWS® Safety Harness



The FitPAWS Safety Harness is a professional-quality canine harness for use with the FitPAWS Peanut and FitPAWS Egg, as well as for other canine rehab use. It is the perfect solution if you are nervous about your dog using the ball or for less-stable canines in rehab. It is also a great tool for helping guide your dog on or off of the Peanut or Egg. The molded rubberized ergonomic handle is lightweight but provides a firm grip so you can comfortably manage your dog's movement, help guide it atop the ball or help guide it down once you've finished training.

The main body is constructed from a top and bottom layer of extremely high-denier, poly-coated fabric which is fully padded inside and triple-stitched at each and every stress point. This offers unrivaled strength along with water resistance. The FitPAWS Safety Harness includes 3 D-Rings conveniently located on the top and both sides. These D-rings can be used for tying leads if the dog is rehabbing on a treadmill.

Our FitPAWS Safety Harness features two adjustable straps of high-strength nylon webbing (underneath the girth and across the chest) for a customized fit. The girth strap also features a safety release. Order correct size by measuring your dog's girth. Use a soft fabric measuring tape and measure your dog around the chest just behind the front legs. Make sure your dog is standing up while you are doing this. You can also circle your dog's chest with a leash and then measure the leash where you mark it.

Product Code

- J1129c
- J1129d
- J1129e
- J1129f
- J1129g

Product/Size/Color

- FitPAWS® Safety Harness/Extra-Large 40" - 42" girth (102cm-107cm)/Red
- FitPAWS® Safety Harness/Large 34" - 40" girth (86cm-102cm)/Yellow
- FitPAWS® Safety Harness/Medium 30" - 34" girth (76cm-86cm)/Blue
- FitPAWS® Safety Harness/Small 26" - 30" girth (66cm-76cm)/Red
- FitPAWS Safety Harness/Extra-Small 18" - 24" girth (46cm-61cm)/Yellow

FPHARN-XL

FPHARN-L

FPHARN-M

FPHARN-S

FPHARN-XS

