

CARPO-FLEX SPORTS WRAPS




designed for use during periods of normal daily activity

USES

- stabilize / support the wrist
- support following removal of a splint or cast
- protect against re-injury
- offload painful arthritic carpal joints
- reduce carpal deviations or overextension

provides light to moderate flexible support for the carpus / wrist



Lines indicate areas that can be trimmed on the wrap to provide a custom fit.

neoprene cuff, elastic straps, and accessory hole for the carpal pad can be trimmed for a custom fit



ITEM

J1004

SIZE

XS - right

J1004 L

XS - left

J1004 A

S - right

J1004 AL

S - left

J1004 B

M - right

J1004 BL

M - left

J1004 C

L - right

J1004 CL

L - left

Carpus-Flex Sports Wraps fasten in a medial direction



XS: 5"W x 3.5" H ; 1.0mm neoprene; 1"W elastic straps (fits small dogs with up to 4.5" wrist circumference)

Dachshund, Fox Terrier, French Bulldog, Jack Russell, King Charles Spaniel, Pug, Sheltie



S: 6"W x 5" H ; 1.5mm neoprene; 1"W and 1.5"W elastic straps (fits small-med dogs with up to 5.5" wrist circumference)

Australian Cattle Dog, Beagle, Border Collie, Brittany, Cocker, Samoyed, Vizsla, Whippet



M: 7"W x 6" H ; 1.5mm neoprene; 1.5"W elastic straps (fits medium dogs with up to 6.5" wrist circumference)

Boxer, Collie, Dalmation, Doberman, Greyhound, Husky, Pointers, Retrievers, Setters, Weimaraner



L: 8"W x 6.75" H ; 2.0mm neoprene; 2"W elastic straps (fits large dogs with up to 7.5" wrist circumference)

German Shepherd, Great Pyrenees, Rottweiler, Great Swiss Mountain Dog, Rhodesian

DORSI-FLEX ASSIST

provides light to moderate assistance for normal hind paw placement

reduces knuckling of the hind paw



Uses:

- reduces knuckling/standing on the dorsum of the paw
- mobility aid for unstable or weak distal hind legs
- assists with tarsal/hock flexion and digital extension
- can be used bilaterally

Features:

- adjustable, elastic tension straps
- includes breathable, anti-skid, Thera-Paw boot
- will not immobilize joints
- lightweight
- easy to put on/take off
- cleans easily by rinsing with warm water

* not to be used for dogs with severe neurological impairments

designed for use as a therapy aid or for daily activity

ITEM #: J1005

SIZE: XS

Includes Thera-Paw size Sm/Med#2 (paw 2.5"W x 3"L) and Tibial cuff up to 4.5" circumference

ITEM #: J1005 A

SIZE: S

Includes Thera-Paw size Med (paw 2.75"W x 3.5"L) and Tibial cuff up to 5.5" circumference

ITEM #: J1005 B

SIZE: M

Includes Thera-Paw size Med/Lrg (paw 3"W x 3.75"L) and Tibial cuff up to 6.5" circumference

ITEM #: J1005 C

SIZE: L

Includes Thera-Paw size Lrg (paw 3.5"W x 4.25"L) Tibial cuff up to 8" circumference