

CONVERSION TABLE

Conversion Table French (Charriere) Gauge to Metric

3	4	5	6	7	8	9	10	11
•	•	•	•	•	•	•	•	•
1mm			2mm			3mm		
12	13	14	15	16	17	18	19	20
•	•	•	•	•	•	•	•	•
4mm			5mm			6mm		
22	24	26	28	30	32	34		
•	•	•	•	•	•	•	•	
	8mm			10mm				
36		38		40		45		
•		•		•		•		
12mm			ACTUAL SIZES			15mm		

THIS CHART IS ACTUAL SIZE