



FitPAWS® Benefits

- Balances out weak areas in the dogs body
 - Improves Core Strength
 - Stabilizes weak areas
 - Improves reaction, control, reflexes and coordination
 - Increases range of motion in joints and elongation of the muscles
 - Improves balance and proprioception
 - Once a dog is on top of FitPAWS Conditioning Equipment, the dogs body will automatically react to the shifting movement, causing the dog to engage different muscle groups simultaneously to remain balanced
 - Modify the balance challenge by increasing or decreasing the inflation inside FitPAWS® products.
-

FitPAWS® Safety Harness

- Keeps dogs safe and in position for controlling movement while rehabbing on top of FitPAWS® products.
- Water resistant Three D-ring Harness is used for rehabbing on land treadmills, underwater treadmills and hydrotherapy pools. (Not a flotation device)
 - 6 sizes for a custom fit: XXS, XS, S, M, L, XL





FitPAWS® Balance Pad

- Low balance and weight bearing tool for dogs recovering from injury or post-operation.
 - Place balance pads side by side, front and back or use 4 FitPAWS Balance Pads for larger dogs
 - Balance pad is made from closed-cell foam, impervious to water for under-water therapy.
-

FitPAWS® Hurdle Set

- Completely adjustable and portable cone and hurdle set includes four 40” bars, 8 cones, 8 height adjustments; from 2”-30”.
- Use for gait training, crawling exercises, cone weaving or jumping.
- Increases coordination, balance, strength, endurance, proprioception, reflexes and gait.
- Encourages use of limbs.
- Includes a black nylon carrying bag.



FitPAWS® Giant Rocker Board

- Use to improve proprioception, weight bearing issues and beginning balance therapy post injury.
- 30X30 padded, non-slip surface
- Easy to clean
- Replacement Mats are available
- Includes insert with exercises



FitPAWS® Peanuts

- Available in 5 sizes: 40cm, 50cm, 60cm, 70cm, & 80cm.
- Improves core strength, increases trunk, stabilizes weak areas, improves balance and proprioception, increases range of motion, and elongation of the muscles.
- Use to assist standing, weight bearing therapy, and rhythmic stabilization.



FitPAWS® Balance Disk & Giant Balance Disk

- Balance and sensory training tool used for core strengthening, weight bearing therapy, and stabilizing weak areas.
- Available in 2 sizes, 14" & 22".
- Smooth side for an increased balance challenge.
- Sensory side provides neurosensory pad stimulation and additional grip ability for grip challenged dogs.
- Use 2 Balance discs side by side to target front and hind end strength. Create an obstacle course using multiple Discs and Balance pads to increase the stability challenge.



FitPAWS® Donut and FitPAWS® Donut Holder

- Lower to the ground; use for weight bearing therapy for large dogs who need to work on front or hind end issues.
 - Multi-Directional movement for an increased stability challenge
 - The FitPAWS® Holder compliments the Donut providing additional stabilization post injury, and to assist beginners.
 - An excellent tool for small dogs!
-

FitPAWS® Paw Pods

- Use for weight bearing therapy on individual limbs, proprioception training, core strengthening and advanced balance training.
- Features sensory bumps for neurosensory stimulation and grip ability.
- Provides two balance challenges: sensory side up is less difficult & challenging, sensory side down is more difficult & challenging.



FitPAWS® Eggs

- Improves balance, coordination and stabilizes weak areas.
 - Use to assist standing, rhythmic stabilization, and balance therapy
 - FitPAWS® Egg offers an increased balance challenge with its' unique shape, smooth surface, and multi-directional movement.
 - Use for advanced proprioception and body awareness.
-

FitBALL® Resistance Bands

Use to add resistance during leash walking, swimming, treadmill work, or sling walking

Can be tied to an individual limb for isolated strength training

Available in 5 resistance levels, XLITE, LITE, MED, HVY & XHVV

Available in precut 4' length, 150 yard bulk rolls, or 12" band loops.

