



## Vet Pro Infusion Pump

“They say great things come in small packages”

The Vet Pro 2000 is specifically designed to be extremely small, but very powerful for today’s veterinary clinic.

User-friendly Numeric Keypad, Full Text LCD

Li-Ion Battery -15 hrs on one charge

Lightweight only 360 gm = < 1 Lb

Adjustable pressure settings


All types of IV solutions even **Blood**

Any type if IV Line can be used  
Performs best with Abbott products



New Piston Pumping Technology for reliability and accuracy +/- 5%

This amazing little pump is lightweight, weighing less than 1 lb and about the size of a soda. The Vet Pro 2000 can operate on A/C power while charging or remove from charging base to operate on battery alone. This will allow for greater portability making this pump the very last pump you will have to or want to buy!



## Vet Pro Twin Infusion Pump

“If you think one is good, two is always better”



The Vet Twin enables you to have the ability to run dual lines at one time on the same pump!

Why tangle with multiple pumps when you can use just ONE!



Vet Twin has all the same specifications as the Vet Pro but is set up for running dual IV sets. This pump will operate off of A/C power while charging battery or on battery for approx 7.5 hrs with a charge time of only 4-6 hours. Extremely small 6 ¾ “ L x 4 ½ “ H x 3” D and lightweight < 1 LB will take care of all clinic needs in one shot!

## Specifications of the Jorvet Vet Pro 2000 Infusion pump

Warranty	1 year limited warranty
Flow rate	0.1 ml/hr to 1200 ml/hr
Types of infusions	Continuous, intermittent, TPN, blood
Accuracy	+/- 5%
KVO Rate	0.1 ml/hr to 5 ml/hr (user adjustable)
Air Sensor	Ultrasonic
Pressure settings	High, Normal, Low- adjustable. Max=10PSI
Power supply	110 – 240 VAC, 50/60 Hz , 20 W
Battery	Li - Ion 7.2 V , 1800 mA Rechargeable
Battery Charging	Automatic when latched in A/C charger
Charge Time	4-6 hours for total charge
Alarms	End Program, Low/End battery, Air inline Door Open, Lock keypad
Mounting	IV Pole mounting or Cage Clip
Dimensions	6 <sup>3</sup> / <sub>4</sub> “ L x 3 <sup>1</sup> / <sub>2</sub> “ W x 3 <sup>3</sup> / <sub>8</sub> “ D
Weight	360 gm (11.6 oz) including mount/charging base
Catalog #	J1097            Vet Pro 2000 J1097M        Vet Pro MRI Compatible J1098            Dual Infusion pump

Made in Israel

Jorgensen Laboratories Inc.  
 1450 N. van Buren Ave  
 Loveland, CO 80538  
 1-800-525-5614  
 Jorvet.com