


Jorvet™



Soft-E-Collar™




The use of Elizabethan collars is widespread in veterinary medicine. There have been a multitude of different designs over the years. Traditional “E” collars have tried to balance the comfort of the patient versus the effectiveness of the “E” collar.

The Soft-E-Collar™ was designed first and foremost for the comfort of the patient, yet it still is effective in offering licking and chewing protection. Clients like the attractive pattern along with the comfortable foam padding.

Pets can move about freely without damaging furniture, and are able to eat and drink more freely than with traditional collars.

The Soft-E-Collar™ is a waterproof heavy nylon material that covers a soft foam doughnut- shaped ring. It is held in place by a strong, flat tie.





Available Sizes

<u>Product Code</u>	<u>Weight (oz.)</u>	<u>Diameter (inches)</u>	<u>Center (Adjustable) (inches)</u>
J-1003A XX-small (cat/dog up to 5 lbs)	2.1	8	1 7/8
J-1003B X-small (cat/dog 5-9 lbs)	2.6	9 1/4	2 3/4
J-1003C Small (cat/dog 10-14 lbs)	3.2	10 3/4	4
J-1003D SS/Medium (cat/dog 15-20 lbs)	5.7	14	4 1/2
J-1003E Small/Medium (dog 20-30 lbs)	7.0	16.5	5
J-1003F Medium (dog 30-55 lbs)	11.3	19	6
J-1003G Large (dog 55-75 lbs)	14.0	23 1/2	7 3/4
J-1003H X-large (dog 75-95 lbs)	19.0	27 1/2	7 1/4
J-1003J XX-large (dog over 95 lbs)	22.00	31 1/2	7 1/4
J-1003S Set of all 9 sizes			

Jorgensen Laboratories, Inc.
 1450 North Van Buren Ave.
 Loveland, CO 80538
 (800) 525-5614
www.jorvet.com
info@jorvet.com

Patent number
 6,044,802